


ICONN


— ORTHOPEDICS —

Rotator Cuff Repair

ICONN Answer II PEEK Suture Anchor
4.75/5.5/6.5

Surgical Protocol by
Geoffrey Connor, MD
Charleston, SC





ICONN's revolutionary approach is to drive value into every aspect of our business; from product design and engineering to distribution and delivery and every step in between.

ICONN Orthopedics
400 Union Hill Drive, Suite 150
Birmingham, AL 35209
800.391.6476

TABLE OF CONTENTS

01	Suture/Anchor/Driver Assembly
03	Arthroscopic Preparation
05	Anchor Pilot Hole Preparation
07	Anchor Insertion
09	Secure Torn Tissue
11	Implant And Instrument List

STEP 1 – SUTURE/ANCHOR/DRIVER ASSEMBLY




FIGURE
01

Anchor/Driver Attachment

Remove ICONN Answer II 4.75MM (147501), 5.5MM (155001), or 6.5MM (165001) Suture Anchor from sterile packaging and place onto ICONN Answer II 4.75/5.5MM Suture Anchor Driver (247506).




FIGURE
02

Suture/Anchor/Driver Attachment

Remove desired number and color combination of ICONN HS Fiber Sutures (HSe-XX-NN) and/or ICONN TetherTape™ (HSXXX) based on surgeon preference. Insert ends of suture(s)/tape through hole in driver shaft near tip (indicated by #1) and pull through to midpoint of suture(s)/tape.




FIGURE
03

Suture/Anchor/Driver Attachment

With driver handle in palm of left hand, **hold trailing ends of suture(s)/tape with left index/thumb**. Using right hand, **apply tension to suture(s)/tape** and wrap around tip of anchor. Join all four strands together and make sure suture(s)/tape are in the suture channel of the anchor.



FIGURE
04

Suture/Anchor/Driver Attachment

Wrap both leading and trailing ends of suture(s)/tape around base of driver handle (indicated by #2 on rear end of driver) so suture(s)/tape are held within slot in driver handle.



FIGURE
05

Excess Suture Wrap

Pull tension on suture(s)/tape and wrap around cleat (under grommet) on top of driver handle (indicated by #'s 3 and 4).

STEP 2 – ARTHROSCOPIC PREPARATION

FIGURE 06

Portal Placement

Standard posterior and anterior portals are utilized for instrument passage and diagnostic arthroscopy. Arthroscope is placed in subacromial space through posterior portal to visually assess pathology.





FIGURE
07

Reduce Cuff

With tissue grasper, reduce cuff to lateral aspect of tuberosity to determine medial-to-lateral mobility of defect.



STEP 3 – ANCHOR PILOT HOLE PREPARATION

FIGURE 08

Select Desired Anchor Position

Place ICONN Answer 4.75/5.5MM Suture Anchor Awl (247502) at a 45° on desired anchor location.





FIGURE
09

Prepare Pilot Hole

Mallet awl to appropriate depth (front of laser line for the 4.75MM anchor and back of the laser line for the 5.5MM and 6.5MM anchors).



STEP 3 – ANCHOR INSERTION

FIGURE 10

Insert Anchor

Remove awl and insert anchor/driver assembly into pilot hold and screw appropriately sized anchor into bone by hand. Continue to advance anchor until the back of the anchor is 1-3mm below the surface of the bone as indicated by the distal end of the laser mark on the driver shaft.





FIGURE
11

Remove Driver

Unwrap suture(s)/tape from cleat on driver handle and pull back on driver to detach from anchor. Confirm anchor is fully seated by direct visualization.



STEP 3 – SECURE TORN TISSUE

FIGURE 12

Secure Sutures

Pass suture(s) through torn tissue using surgeon's preferred method. Secure torn tissue using standard arthroscopic knot tying techniques.





FIGURE 13


Repeat

Repeat previous steps until desired number of anchors are inserted and torn tissue is fully secured. Upon completion close arthroscopic incisions.




IMPLANT AND INSTRUMENT LIST

PART NUMBER	DESCRIPTION
147501	ICONN Answer II 4.75MM Suture Anchor
155001	ICONN Answer II 5.5MM Suture Anchor
165001	ICONN Answer II 6.5MM Suture Anchor
200001 & 200002	ICONN Answer Instrument Tray Base and Lid
247506	ICONN Answer II 4.75/5.5MM Suture Anchor Driver
247502	ICONN Answer 4.75/5.5MM Suture Anchor Awl
HSe-2C-NN	HS Fiber White/Blue #2 Suture
HSe-2K-NN	HS Fiber White/Black #2 Suture
HSe-2G-NN	HS Fiber White/Green #2 Suture
HSe-2XC-NN	HS Fiber Blue #2 Suture
HS226	1.6mm TetherTape™ White/Blue
HS227	1.6mm TetherTape™ White/Black
HS228	1.6mm TetherTape™ White/Blue/Blue/Black
HS229	1.6mm TetherTape™ White/Green



We are designing products, and even more so,
a company which is equally as innovative from
a performance standpoint as from a
cost-reduction standpoint.



ICONN Orthopedics
400 Union Hill Drive, Suite 150
Birmingham, AL 35209
800.391.6476