


# ICONN

## ORTHOPEDICS

### ICONN Answer PEEK Suture Anchor Product Information



#### ICONN Answer 4.75 & 5.5MM Suture Anchors

- Rotator Cuff Repair
- Loaded with 2 or 3 #2 HS Fiber Sutures
- 5.5MM has identical thread geometry to 4.75MM to act as a “rescue” anchor in soft bone


#### HS Fiber™

- Ultra High Molecular Weight Polyethylene for maximum strength
- 4 color options to aid in suture management
- 10x stronger than steel on a weighted basis


- ✓ Designed from the ground up with a dual focus on quality and cost, resulting in implants priced at 1/2 to 1/3 of competitive anchors
- ✓ Fully threaded “corkscrew”-tapered design for maximum bone purchase
- ✓ Wrap-around suture channel eliminates eyelet failure and allows for ease of assembly

- ✓ Reusable driver eliminates waste and provides a higher quality insertion tool
- ✓ Implants/sutures individually packaged sterile, reducing inventory cost
- ✓ In-surgery assembly for more efficient delivery

#### Pullout Strength After Cyclic Loading



#### HS Fiber™ Strength Comparison







ICONN's revolutionary approach is to drive value  
into every aspect of our business; from product  
design and engineering to distribution and delivery  
and every step in between. We are not simply selling  
a lower cost product. We are designing products,  
and even more so, a company which is equally as  
innovative from a performance standpoint as from a  
cost-reduction standpoint.

ICONN Orthopedics  
400 Union Hill Drive, Suite 150  
Birmingham, AL 35209  
800.391.6476